## 40" Ductless Fume Hood

## Model # SS-340-DCH

The **40" Wide Ductless Fume Hood** combines powerful air flow with a sturdy and compact enclosure for optimal respiratory and environmental protection. This fume hood helps protect the operator and environment from hazardous fumes or particulate generated from applications performed within the hood.

The **40" Ductless Fume Hood's** powerful fan directs fumes into the filter chamber containing high quality filter media selected depending on the application. Filter media choices include **HEPA** [up to 99.97% efficiency on particles 0.3 microns and larger], **ASHRAE** [up to 95% efficiency on particles 0.5 microns and larger, **Activated Carbon**, and **specialty-blended filtration** [i.e. acid gas, mercury, aldehyde, and ammonia].

Typical uses for the **40" Ductless Fume Hood** include chemical fumes, pharmaceutical compounding, soldering applications, dust removal, biological applications, and a variety of applications requiring fume and dust removal. The **40" Ductless Fume Hood** can provide an effective and economical solution for many commercial, medical, laboratory, and industrial applications. For additional visibility, a transparent version is available at extra cost.

Our line of ductless fume hoods is not intended for sterile compounding applications.



MADE IN

THE USA

40" Ductless Fume Hood



40" DCH provides fume control for a variety of medical and lab applications.

## **Product Features**

- Includes: variable speed controller (115V only), velometer, and LED light
- **Optional:** clear anti-static vinyl curtains, Magnehelic Gauge, hour counter, and spill tray
- **Standard Sizes:** 12", 18", 24", 30", 40", 50", 60", and 70" (custom sizes available)
- · Quiet operation
- · Long filter life
- · Simple, quick "no tool" filter change
- · Portable and lightweight
- Reliable, low maintenance operation



40" DCH provides chemical fume control for use with various substances.



21221 FM 529 Rd, Cypress, TX 77433 1.800.799.4609 • sales@sentryair.com

## **Product Specifications**





Shown with outer dimensions. Dimensions are approximate.



21221 FM 529 Rd Cypress, TX 77433 1.800.799.4609 sales@sentryair.com

www.sentryair.com

Hood Dimensions	Outer: 24" x 42" x 25"
(D x W x H)	<b>Interior:</b> 23.5" x 40" x 24"
(DAWAII)	Access Area: 8" x 40"
Hood Material	Sides: 0.25" white acrylic
	Base: 0.5" white high
	density polyethylene
	<b>Top:</b> 0.5" clear
	polycarbonate
	Front Window: 0.25"
Bloom Bloom Store	clear polycarbonate
<b>Blower Dimensions</b>	15.5"L x 13" W x 14.5" H
	16 Ga. carbon steel with
Blower Material	polyester epoxy coating
	on all surfaces
Weight	115 lbs
Air Volume	350 CFM
(approximately)	(varies with filter media)
	Carbon Pre-filter + HEPA
	Filter – 130 FPM high,
Inlet Velocities on	25 FPM low
High (approximately)	Carbon Pre-filter + Heavy
	Duty Carbon Filter –
	115 FPM high, 20 FPM low
Electrical	115/1/60, 2 amps
	220/1/50, 0.75 amps
	8' grounded power cord
	with NEMA 5-15P plug
Sound Level	61 to 66 dba (low to high)
Maximum	Operating temperature
Temperature	not to exceed 170° F
	<b>HEPA</b> (up to 99.97%
Main Filter Options	efficient on particles
	down to 0.3 microns)
	ASHRAE (up to 95%
	efficient on particles
	down to 0.5 microns)
	Activated Carbon
	Specialty-Blended Filter
	Media (i.e. acid gas,
	mercury, aldehyde, and
	ammonia)
Warranty	Limited two-year warranty
	from date of shipment on
	defects due to materials
	or workmanship.
Patent #	5,843,197